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Cool Coil and Fans

Hot Coil and Fans, Compressor & Condenser

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
As can be seen from the experiment, a normal probe, due to its vulnerability to intermittent air currents, shows high temperature fluctuations, especially when there are door openings or during defrosts cycles. With the standard deviation and the graphs in the experiment, we can see without a doubt that temperature measurements recorded by a normal probe are very unstable. On the other hand, a glycol-encased probe consistently shows stable temperature measurements and undoubtedly is better at representing the temperature of vaccines inside refrigerators. With the experiment being tested with different types of refrigerators, this experiment not only meticulously demonstrated the very difference between glycol-encased probes and normal probes, but also corroborated the CDC's recommendation. Using a normal probe can result in many problems, of which the most obvious is wasting viable vaccines due to inaccurate readings. Therefore, choosing a glycol-encased probe over a normal probe would not only provide more accurate temperature readings and help to protect your medical materials, but also save you a lot of money in the long term.

## References

Centers for Disease Control and Prevention. ["Vaccine Storage & Handling Toolkit."](#) June 2016. 15 July 2016

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| Label | Name                                            | Range                                | °C/°F |     |
|-------|-------------------------------------------------|--------------------------------------|-------|-----|
| Set   | Set point                                       | LS÷US                                | 3.0   | --- |
| Hy    | Differential                                    | 0.1÷25.5°C/ 1÷ 255°F                 | 2.0   | Pr1 |
| LS    | Minimum set point                               | -50°C÷SET/-58°F÷SET                  | -50.0 | Pr2 |
| US    | Maximum set point                               | SET÷110°C/ SET ÷ 230°F               | 110   | Pr2 |
| Ot    | Thermostat probe calibration                    | -12÷12°C /-120÷120°F                 | 0.0   | Pr1 |
| P3P1  | Third probe presence                            | n=not present; Y=pres.               | n     | Pr2 |
| O31   | Third probe calibration                         | -12÷12°C /-120÷120°F                 | 0     | Pr2 |
| P4P   | Fourth probe presence                           | n=not present; Y=pres.               | n     | Pr2 |
| O4    | Fourth probe calibration                        | -12÷12°C /-120÷120°F                 | 0     | Pr2 |
| OdS   | Outputs delay at start up                       | 0÷255 min                            | 0     | Pr2 |
| AC    | Anti-short cycle delay                          | 0 ÷ 50 min                           | 1     | Pr1 |
| CCt   | Continuous cycle duration                       | 0.0÷24.0h                            | 0.0   | Pr2 |
| CCS   | Set point for continuous cycle                  | (-55.0÷150.0°C) (-67÷302°F)          | 3     | Pr2 |
| COu   | Compressor ON time with faulty probe            | 0 ÷ 255 min                          | 15    | Pr2 |
| COF   | Compressor OFF time with faulty probe           | 0 ÷ 255 min                          | 30    | Pr2 |
| CH    | Kind of action                                  | CL=cooling; Ht= heating              | cL    | Pr1 |
| CF    | Temperature measurement unit                    | °C ÷ °F                              | °C    | Pr2 |
| rES   | Resolution                                      | in=integer; dE= dec.point            | dE    | Pr1 |
| dLy   | Display temperature delay                       | 0 ÷ 20.0 min (10 sec.)               | 0     | Pr2 |
| IdF   | Interval between defrost cycles                 | 1 + 120 ore                          | 8     | Pr1 |
| MdF   | (Maximum) length for defrost                    | 0 + 255 min                          | 20    | Pr1 |
| dFd   | Displaying during defrost                       | rt, it, SET, DEF                     | it    | Pr2 |
| dAd   | MAX display delay after defrost                 | 0 ÷ 255 min                          | 30    | Pr2 |
| ALc   | Temperat. alarms configuration                  | rE= related to set;<br>Ab = absolute | Ab    | Pr2 |
| ALU   | MAXIMUM temperature alarm                       | Set÷110.0°C; Set÷230°F               | 110   | Pr1 |
| ALL   | Minimum temperature alarm                       | -50.0°C÷Set/ -58°F÷Set               | -50.0 | Pr1 |
| AFH   | Differential for temperat. alarm recovery       | (0.1°C÷25.5°C) (1°F÷45°F)            | 1     | Pr2 |
| ALd   | Temperature alarm delay                         | 0 ÷ 255 min                          | 15    | Pr2 |
| dAo   | Delay of temperature alarm at start up          | 0 ÷ 23h e 50'                        | 1.3   | Pr2 |
| AL2   | Condenser for low temperat. alarm               | (-55 ÷ 150°C) (-67÷ 302°F)           | -40   | Pr2 |
| AU2   | Condenser for high temperat. alarm              | (-55 ÷ 150°C) (-67÷ 302°F)           | 110   | Pr2 |
| AH2   | Differ. for condenser temp. alar. recovery      | [0.1°C ÷ 25.5°C] [1°F ÷ 45°F]        | 5     | Pr2 |
| Ad2   | Condenser temperature alarm delay               | 0 ÷ 254 (min.) , 255=nU              | 15    | Pr2 |
| dA2   | Delay of cond. temper. alarm at start up        | 0.0 ÷ 23h 50'                        | 1.3   | Pr2 |
| bLL   | Compr. off for condenser low temperature alarm  | n(0) - Y(1)                          | n     | Pr2 |
| AC2   | Compr. off for condenser high temperature alarm | n(0) - Y(1)                          | n     | Pr2 |
| i1P   | Digital input polarity                          | oP=opening;CL=closing                | cL    | Pr1 |
| i1F   | Digital input configuration                     | EAL, bAL, PAL, dor, dEF, Htr, AUS    | dor   | Pr1 |
| diid  | Digital input alarm delay                       | 0÷255min                             | 15    | Pr1 |
| nPS   | Number of activation of pressure switch         | 0 ÷ 15                               | 15    | Pr2 |
| odc   | Compress status when open door                  | no; Fan; CPr; F_C                    | no    | Pr2 |
| rrd   | Regulation restart with door open alarm         | n - Y                                | y     | Pr2 |
| HES   | Differential for Energy Saving                  | (-30°C÷30°C) (-54°F÷54°F)            | 0     | Pr2 |
| Adr   | Serial address                                  | 0÷247                                | 1     | Pr2 |
| PbC   | Kind of probe                                   | Ptc; ntc                             | ntc   | Pr1 |
| onF   | on/off key enabling                             | nu, oFF; ES                          | nu    | Pr2 |
| dP1   | Room probe display                              | --                                   | --    | Pr1 |
| dP3   | Third probe display                             | --                                   | --    | Pr1 |
| dP4   | Fourth probe display                            | --                                   | --    | Pr2 |
| rSE   | Real set point value                            | actual set                           | --    | Pr2 |
| rEL   | Software release                                | --                                   | --    | Pr2 |
| Ptb   | Map code                                        | --                                   | --    | Pr2 |

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